

CERTIFICATE OF COMPLIANCE



*TNO Certification
hereby declares that the product*

Flow Sensor type T series, size range T025, T050, T075, T100 and T150

From

Micro Motion, Boulder (Co), USA

*has been evaluated for compliance with the Hygienic Equipment Design Criteria
of the EHEDG, Document No. 8, by:*

*TNO Nutrition and Food Research at Zeist, Netherlands
and meets the criteria of this document as demonstrated by:*

Evaluation Report No. V99.663

Signed  _____
Evaluation Officer

Date January 10 2002

Signed  _____
Managing Director, TNO Certification

Date January 10 2002



Certificate No. COI-3682

*TNO Certification BV, P.O. Box 541, 7300 AM Apeldoorn, Netherlands
©EHEDG*