

The Industry's Most Complete Line of Precision Cleaning Equipment

Ultrasonic Components



Ultrasonic Tanks

The cornerstone of any ultrasonic system, Branson makes standard stainless steel ultrasonic tanks ranging from 5 to 100 gallons, and custom tanks in nearly any size. A number of options are available, including filter/recirculation packages and parts baskets.

Generators

Available in 3 model types and 5 frequencies (25, 40, 80, 120 and 170 kHz) offering a variety of features to match your cleaning needs.

Immersible Transducers

Branson 316L stainless steel immersible transducers provide versatile ways to add ultrasonics to new or existing tanks. They are available in standard and custom sizes. Also available in Hastelloy[†].

Support Components

Hot Air Dryers

Temperature-controlled heaters and recirculating hot air provide up to 250°F drying. Options include high-performance blowers, power operated cover, air-blow-off nozzles and HEPA filtration.

Coalescer

Designed to separate oily soils from aqueous cleaning solutions, extending the solution's useful life.

Aquafier

This closed-loop industrial water purification system produces high-quality deionized (DI) rinse water at a rate of up to 5 gpm.

